





NEUROSAFE

Glutton™ Aspiration Catheter

The Glutton™ Aspiration Catheter is intended for use in the revascularization of patients with acute ischemic stroke secondary to intracranial large vessel occlusive disease (with the internal carotid, middle cerebral – M1 and M2 segments, basilar, and vertebral arteries) within 8 hours of symptom onset.



LARGE LUMEN **0.072"**
for aspiration, capture large clot



Large Lumen
Ultra Flexibility
Great Efficiency

Hybrid Structure of Coil and Braiding,
BALANCED KINK RESISTANCE AND TRACKABILITY.

FEATURES & BENEFITS

- - Hybrid structure of coil and braiding, combined flexibility and maneuverability;
- - Larger inner diameter lumen, higher aspirating efficiency;
- - Distal multi-segment design of gradual hardness;
- - Proximal composite polymer enhanced pushability;
- - Compatible with 6F long sheath;
- - Hydrophilic coating on the distal segment easily reach distal vessels / occlusions.

Order Information:

| Catalog No. | Fr. Size | Inner Diameter (inch) | Outer Diameter (inch) | Effective Length (cm) |
|-------------|----------|-----------------------|-----------------------|-----------------------|
| 5F-115 | 5F | 0.058 | 0.070 | 115 |
| 5F-125 | 5F | 0.058 | 0.070 | 125 |
| 5F-132 | 5F | 0.058 | 0.070 | 132 |
| 6F-105 | 6F | 0.072 | 0.084 | 105 |
| 6F-115 | 6F | 0.072 | 0.084 | 115 |
| 6F-125 | 6F | 0.072 | 0.084 | 125 |
| 6F-115S | 6F | 0.072 | 0.084 | 115 |
| 6F-125S | 6F | 0.072 | 0.084 | 125 |
| 6F-132S | 6F | 0.072 | 0.084 | 132 |

NeuroSafe Medical Co., Ltd.

Address: Building B, No. 10, Keji 1st Rd, Hi-Tech Industrial Development Zone, Zhuhai, Guangdong, China

Tel: +86 (756) 3635 345 | Email: sales@neurosafe.cn | Post code: 519 000

Website: www.neurosafemed.com